

Product Datasheet

Hemoglobin (Hb) Antibody HM061

Product Name	Hemoglobin (Hb) Antibody
Catalog Number	HM061
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified Human hemoglobin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Hemoglobin
Matched Pairs	HM754 or HM753 (capture)
Preservatives	0.09% Sodium Azide