

Product Datasheet

Hemoglobin (Hb) Antibody HM061

Product Name	Hemoglobin (Hb) Antibody
Catalog Number	HM061
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA, LF
Immunogen	Purified human hemoglobin
Brand	EastCoast Bio
Storage	Short Term: 2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	≥ 95% by SDS-PAGE
Buffer	Phosphate buffered saline, pH 7.0 ± 0.2 with 0.09% sodium azide as a preservative.
Specificity	Hemoglobin (Hb)
Matched Pairs	HM754 or HM753 (capture)
Preservatives	0.09% Sodium Azide