

Product Datasheet

Transferrin Antibody HM063

Product Name	Transferrin Antibody
Catalog Number	HM063
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified human transferrin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Matched Pairs	HM065 (detection)
Preservatives	0.09% Sodium Azide