

Product Datasheet

NT-proBNP (13-27) Antibody HM071

Product Name	NT-proBNP (13-27) Antibody
Catalog Number	HM071
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified recombinant human NT-proBNP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	NT-proBNP
Matched Pairs	HM074 or HM316 (detection) and LA313
Preservatives	0.09% Sodium Azide