

Product Datasheet

Procalcitonin (PCT-calcitonin) Antibody HM081

Product Name	Procalcitonin (PCT-calcitonin) Antibody
Catalog Number	HM081
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified recombinant human Procalcitonin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	PCT-calcitonin epitope
Matched Pairs	HM078 (capture)
Preservatives	0.09% Sodium Azide