

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody HM082

Product Name	Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody
Catalog Number	HM082
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Native HIV antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	HIV-1 p24
Matched Pairs	JP040 or HM1202 (detection)