

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody HM083

Product Name	Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody
Catalog Number	HM083
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Native HIV antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20 °C
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline pH 7.4 ± 0.2, with 0.15M NaCl
Specificity	HIV-1 p24
Matched Pairs	HM084 or JP040
Preservatives	0.1 % Proclin 300