

Product Datasheet

Malaria (pf-HRP-2) Antibody HM087

Product Name	Malaria (pf-HRP-2) Antibody
Catalog Number	HM087
Species/Host	Mouse
Tested Applications	FC, ELISA, Lateral Flow
Immunogen	Recombinant Malaria HRP-2 antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Plasmodium falciparum (malaria) HRP-2
Matched Pairs	HM086 (capture)
Preservatives	0.1% Proclin 300