

Product Datasheet

Malaria (pf-HRP-2) Antibody HM087

Product Name	Malaria (pf-HRP-2) Antibody
Catalog Number	HM087
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, FC
Immunogen	Recombinant Malaria HRP-2 antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Plasmodium falciparum (malaria) HRP-2
Matched Pairs	HM086 (capture)
Preservatives	0.1% Proclin 300