

Product Datasheet

Human Immunoglobulin M (IgM) Antibody HM098

Product Name	Human Immunoglobulin M (IgM) Antibody
Catalog Number	HM098
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Human IgM
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Human IgM
Preservatives	0.1% Proclin 300