

## Product Datasheet

### Pregnancy Associated Plasma Protein A (PAPP-A) Antibody HM1015

<b>Product Name</b>	Pregnancy Associated Plasma Protein A (PAPP-A) Antibody
<b>Catalog Number</b>	HM1015
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Purified PAPP-A
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95% by SDS-PAGE
<b>Associated Products</b>	Native PAPP-A antigen 496-90
<b>Buffer</b>	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
<b>Specificity</b>	PAPP-A
<b>Matched Pairs</b>	HM1014 (capture)
<b>Preservatives</b>	0.05% Sodium Azide