

Product Datasheet

Pregnancy Associated Plasma Protein A (PAPP-A) Antibody HM1015

Product Name	Pregnancy Associated Plasma Protein A (PAPP-A) Antibody
Catalog Number	HM1015
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Purified PAPP-A
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Associated Products	Native PAPP-A antigen 496-90
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	PAPP-A
Matched Pairs	HM1014 (capture)
Preservatives	0.05% Sodium Azide