

Product Datasheet

Estrone-3-Sulfate (E1S) Antibody HM1016

Product Name	Estrone-3-Sulfate (E1S) Antibody
Catalog Number	HM1016
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Estrone-3-Sulfate-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Cross Reactivity	Estriol <6%; Estriol-3-glucuronide, 17 β -Estradiol, Estradiol-3-glucuronide <5%; 17 α -Ethinylestradiol <3%; Estrone-3-glucuronide, 17 α -hydroxyprogesterone, Progesterone, Pregnenediol-3-glucuronide, Testosterone, Aldosterone, Androstenedione <1%
Specificity	Estrone and Estrone-3-sulfate
Matched Pairs	LA550 or LA551 (capture)
Preservatives	0.05% Sodium Azide