

Product Datasheet

Wesselsbron Virus (WSLV-NS1) Antibody HM1017

Product Name	Wesselsbron Virus (WSLV-NS1) Antibody
Catalog Number	HM1017
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant Wesselsbron Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1023 (capture or detection)
Preservatives	0.05% Sodium Azide