

Product Datasheet

Wesselsbron Virus (WSLV-NS1) Antibody HM1023

Product Name	Wesselsbron Virus (WSLV-NS1) Antibody
Catalog Number	HM1023
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Wesselsbron Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1017 (capture or detection)
Preservatives	0.05% Sodium Azide