

## Product Datasheet

### Insulin-Like Growth Factor Binding Protein 1 (IGFBP-1) Antibody HM1025

<b>Product Name</b>	Insulin-Like Growth Factor Binding Protein 1 (IGFBP-1) Antibody
<b>Catalog Number</b>	HM1025
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Purified Human IGFBP-1
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0
<b>Specificity</b>	IGFBP-1
<b>Matched Pairs</b>	HM1024 (capture)
<b>Preservatives</b>	0.09% Sodium Azide