

Product Datasheet

Insulin-Like Growth Factor Binding Protein 1 (IGFBP-1) Antibody HM1025

Product Name	Insulin-Like Growth Factor Binding Protein 1 (IGFBP-1) Antibody
Catalog Number	HM1025
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Human IGFBP-1
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	IGFBP-1
Matched Pairs	HM1024 (capture)
Preservatives	0.09% Sodium Azide