

Product Datasheet

African Swine Fever Virus (ASFV-p30) Antibody HM1031

Product Name	African Swine Fever Virus (ASFV-p30) Antibody
Catalog Number	HM1031
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant ASFV strain Morara/Georgia
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1029 or HM1027 (detection)
Preservatives	0.05% Sodium Azide