

Product Datasheet

Ebola Virus (Zaire-VP40) Antibody HM1032

Product Name	Ebola Virus (Zaire-VP40) Antibody
Catalog Number	HM1032
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Recombinant Ebola Virus Zaire
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Ebola Virus Zaire
Preservatives	0.05% Sodium Azide