

## Product Datasheet

### Pancreatic Elastase Antibody HM1033

<b>Product Name</b>	Pancreatic Elastase Antibody
<b>Catalog Number</b>	HM1033
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	+2-8 °C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95 % by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2 ± 0.2
<b>Cross Reactivity</b>	No reactivity with pig elastase; pig, bovine, human hemaglobin; pig, bovine, human transferrin; human calprotectin; bovine, human lactoferrin
<b>Specificity</b>	Pancreatic Elastase
<b>Matched Pairs</b>	HM1034 (capture)
<b>Preservatives</b>	0.09% Sodium Azide