

## Product Datasheet

### Pancreatic Elastase Antibody HM1033

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Pancreatic Elastase Antibody  |
| <b>Catalog Number</b>      | HM1033  |
| <b>Species/Host</b>        | Mouse   |
| <b>Tested Applications</b> | ELISA   |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody  |
| <b>Storage</b>             | 2-8°C   |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Purity</b>              | >95% by SDS-PAGE  |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.2   |
| <b>Cross Reactivity</b>    | No reactivity with pig elastase; pig, bovine, human hemaglobin; pig, bovine, human transferrin; human calprotectin; bovine, human lactoferrin |
| <b>Specificity</b>         | Pancreatic Elastase   |
| <b>Matched Pairs</b>       | HM1034 (capture)  |
| <b>Preservatives</b>       | 0.09% Sodium Azide  |