

Product Datasheet

Hydrocodone Antibody HM1036

Product Name	Hydrocodone Antibody
Catalog Number	HM1036
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Hydrocodone-BSA
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Morphine, Codeine, Hydromorphone, 6-MAM, Oxymorphone <1%
Specificity	Hydrocodone & Oxycodone
Matched Pairs	LA289 (capture)
Preservatives	0.05% Sodium Azide