

# Product Datasheet

## Human CD126 ELISA Kit 850.510

Product Name	Human CD126 ELISA Kit
Catalog Number	850.510
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	5 pg/ml
Tested Applications	ELISA
Alternative Names	gp80 /IL-6R
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble receptors. No interference with neither soluble gp130 nor IL-6.
Specificity	Recognizes both natural and recombinant human soluble IL-6R
Sample Types	Serum, Plasma, Cell culture supernatant