

Product Datasheet

Lysergic acid diethylamide (LSD) Antibody HM1045

Product Name	Lysergic acid diethylamide (LSD) Antibody
Catalog Number	HM1045
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	LSD-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	0.015M Potassium phosphate, pH 7.2 with 0.85% NaCl
Specificity	LSD
Preservatives	0.05% Sodium Azide