

Product Datasheet

Pancreatic Lipase Antibody, Canine HM1046

Product Name	Pancreatic Lipase Antibody, Canine
Catalog Number	HM1046
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Immunogen	Native canine pancreatic lipase
Brand	EastCoast Bio
Form/Appearance	Purified monoclonal antibody
Storage	Short Term: 2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	Mouse IgG
Clonality	Monoclonal
Purity	≥ 95 % by SDS-PAGE
Buffer	Phosphate-buffered saline, pH 7.0 ± 0.2, 0.02% sodium azide as a preservative.
Matched Pairs	HM1047 (detection)
Preservatives	0.02% Sodium Azide