

Product Datasheet

Pancreatic Lipase Antibody, Canine HM1047

Product Name	Pancreatic Lipase Antibody, Canine
Catalog Number	HM1047
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Native Canine Pancreatic Lipase
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Matched Pairs	HM1046 (capture)
Preservatives	0.02% Sodium Azide