

Product Datasheet

Human CD130 ELISA Kit 850.520.048

Product Name	Human CD130 ELISA Kit
Catalog Number	850.520.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	56.25-1800 pg/ml
Sensitivity	15 pg/ml
Tested Applications	ELISA
Alternative Names	gp130
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble receptors. No interference with neither soluble gp80 nor IL-6.
Specificity	Recognizes both natural and recombinant human soluble gp130
Sample Types	Serum, Plasma, Cell culture supernatant