

Product Datasheet

Methylone Antibody HM1067

Product Name	Methylone Antibody
Catalog Number	HM1067
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Methylone-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	4-Chloromethcathinone HCl 100%; 4-Bromomethcathinone HCl 90%; Mephedrone 50%; 4-Methylephedrine HCl, 4-Methylethcathinone 5%; Methamphetamine 3%; MDMA, Amphetamine, MDPV, Pentedrone, Pentylone, a-PVP <1%
Specificity	Methylone
Matched Pairs	LA608 (capture)
Preservatives	0.05% Sodium Azide