

Product Datasheet

Human Immunoglobulin G (IgG) Antibody HM1070

Product Name	Human Immunoglobulin G (IgG) Antibody
Catalog Number	HM1070
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Human IgG
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Human IgA, IgM, IgE
Specificity	Human IgG (IgG1, IgG2a, IgG3, IgG4)
Preservatives	0.09% Sodium Azide