

# Product Datasheet

## Red Blood Cell (RBC), Human Antibody HM1079

Product Name	Red Blood Cell (RBC), Human Antibody
Catalog Number	HM1079
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Human RBC
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term ( $\leq$ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Preservatives	0.05% Sodium Azide