

Product Datasheet

African Swine Fever Virus (ASFV-pp62) Antibody HM1087

Product Name	African Swine Fever Virus (ASFV-pp62) Antibody
Catalog Number	HM1087
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant ASFV strain BA71V
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	ASFV pp62 protein
Matched Pairs	HM1080 or HM1086 (detection)
Preservatives	0.05% Sodium Azide