

Product Datasheet

Neurofilament L (NfL) Antibody HM1095

Product Name	Neurofilament L (NfL) Antibody
Catalog Number	HM1095
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant NfL protein (insect cell expressed)
Alternative Names	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS). Amyotrophic lateral sclerosis (ALS)
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Pair with HM1096 or HM1184 for best results in Lateral Flow
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1096, HM1098, HM1101 or HM1184 (detection)
Preservatives	0.05% Sodium Azide