

## Product Datasheet

### Neurofilament L (NfL) Antibody HM1095

<b>Product Name</b>	Neurofilament L (NfL) Antibody
<b>Catalog Number</b>	HM1095
<b>Description</b>	Ion Exchange Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	Recombinant NfL protein (insect cell expressed)
<b>Alternative Names</b>	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS). Amyotrophic lateral sclerosis (ALS)
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	Pair with HM1096 or HM1184 for best results in Lateral Flow
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
<b>Matched Pairs</b>	HM1096, HM1098, HM1101 or HM1184 (detection)
<b>Preservatives</b>	0.05% Sodium Azide