

Product Datasheet

Neurofilament L (NfL) Antibody HM1098

Product Name	Neurofilament L (NfL) Antibody
Catalog Number	HM1098
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant NfL protein (insect cell expressed)
Alternative Names	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS). Amyotrophic lateral sclerosis (ALS)
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1095 (capture)
Preservatives	0.05% Sodium Azide