

Product Datasheet

Neurofilament L (NfL) Antibody
HM1101

| | |
|---------------------|--|
| Product Name | Neurofilament L (NfL) Antibody |
| Catalog Number | HM1101 |
| Description | Ion Exchange Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Recombinant NfL protein (insect cell expressed) |
| Alternative Names | Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS). Amyotrophic lateral sclerosis (ALS) |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Note | Pair with HM1095 for best results in Lateral Flow |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Matched Pairs | HM1095, HM1097 or HM1100 (detection) |
| Preservatives | 0.05% Sodium Azide |