

Product Datasheet

Neurofilament L (NfL) Antibody HM1101

Product Name	Neurofilament L (NfL) Antibody
Catalog Number	HM1101
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant NfL protein (insect cell expressed)
Alternative Names	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS). Amyotrophic lateral sclerosis (ALS)
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Pair with HM1095 for best results in Lateral Flow
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1095, HM1097 or HM1100 (detection)
Preservatives	0.05% Sodium Azide