

Product Datasheet

Helicobacter pylori Antibody HM1107

Product Name	Helicobacter pylori Antibody
Catalog Number	HM1107
Species/Host	Mouse
Tested Applications	ELISA, IFA
Immunogen	Purified Helicobacter pylori
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Cross Reactivity	Not reactive with Campylobacter, Cryptosporidium, E. coli, Enterococcus, Giardia, Salmonella, Shigella
Specificity	H. pylori flagellin
Matched Pairs	J353 (capture)
Preservatives	0.1% Sodium Azide