

Product Datasheet

Clostridium difficile Glutamate Dehydrogenase (GDH) Antibody HM1110

Product Name	Clostridium difficile Glutamate Dehydrogenase (GDH) Antibody
Catalog Number	HM1110
Species/Host	Mouse
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	C. difficile GDH
Matched Pairs	HM1108 or HM1109 (capture)
Preservatives	0.1% Sodium Azide