

Product Datasheet

Human CD102 ELISA Kit 850.550.048

Product Name	Human CD102 ELISA Kit
Catalog Number	850.550.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.75-24 U/ml
Sensitivity	0.2 U/ml
Tested Applications	ELISA
Alternative Names	ICAM-2
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human adhesion molecules
Specificity	Recognizes both natural and recombinant human soluble ICAM-2
Sample Types	Serum, Plasma, Cell culture supernatant