

Product Datasheet

Creatine Kinase B (CK-B) Antibody HM1128

Product Name	Creatine Kinase B (CK-B) Antibody
Catalog Number	HM1128
Description	Creatine Kinase B (CK-B) Antibody
Species/Host	Mouse
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	>98%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM144 (capture)
Preservatives	0.05% Sodium Azide