

Product Datasheet

Creatine Kinase MB (CK-MB) Antibody HM1128

Product Name	Creatine Kinase MB (CK-MB) Antibody
Catalog Number	HM1128
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	≥ 98% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 ± 0.2
Matched Pairs	HM144 (capture)
Preservatives	0.05% Sodium Azide