

Product Datasheet

Fentanyl Antibody HM1132

Product Name	Fentanyl Antibody
Catalog Number	HM1132
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Fentanyl-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Acetyl Norfentanyl 90%; Acetyl Fentanyl 60%; Norfentanyl & Alfentanil 30%; Carfentanil 1%
Specificity	Fentanyl
Matched Pairs	LA395 (capture)
Preservatives	0.05% Sodium Azide