

Product Datasheet

Pepsin Antibody HM1134

Product Name	Pepsin Antibody
Catalog Number	HM1134
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Human Pepsin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0 ± 0.2
Matched Pairs	HM1136 (capture)
Preservatives	0.09% Sodium Azide