

# Product Datasheet

## Human CD50 ELISA Kit 850.560.096

Product Name	Human CD50 ELISA Kit
Catalog Number	850.560.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.625-20 ng/ml
Sensitivity	0.3 ng/ml
Tested Applications	ELISA
Alternative Names	ICAM-3
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human adhesion molecules
Specificity	Recognizes both natural and recombinant human soluble ICAM-3
Sample Types	Serum, Plasma, Cell culture supernatant