

Product Datasheet

Human CD44 ELISA Kit 850.570.048

Product Name	Human CD44 ELISA Kit
Catalog Number	850.570.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.31-10 ng/ml
Sensitivity	0.113 ng/ml
Tested Applications	ELISA
Alternative Names	PGP-1, H-CAM
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human adhesion molecules
Specificity	Recognizes both natural and recombinant human soluble Pgp-1
Sample Types	Serum, Plasma, Cell culture supernatant