

Product Datasheet

Risperidone Antibody HM1152

Product Name	Risperidone Antibody
Catalog Number	HM1152
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Risperidone-BSA
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Amitriptyline, Aripiprazole, Clomipramine, Clonapine, Cyclobenzaprine, Desipramine, Doxepin, Imipramine, Maprotiline, Nordoxepin, Nortriptyline, Olanzapine, Quetiapine, Trimipramine, Ziprasidone <1%
Specificity	Risperidone
Matched Pairs	LA656 (capture)
Preservatives	0.05% Sodium Azide