

Product Datasheet

Severe Fever with Thrombocytopenia Syndrome Virus (SFTS-NP) Antibody HM1155

| | |
|---------------------|---|
| Product Name | Severe Fever with Thrombocytopenia Syndrome Virus (SFTS-NP) Antibody |
| Catalog Number | HM1155 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Recombinant SFTS Virus nucleoprotein |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | >90% by SDS-PAGE |
| Buffer | 0.15M Potassium Phosphate, 0.85% NaCl, pH 7.2 |
| Preservatives | 0.05% Sodium Azide |