

## Product Datasheet

### Severe Fever with Thrombocytopenia Syndrome Virus (SFTS-NP) Antibody HM1156

<b>Product Name</b>	Severe Fever with Thrombocytopenia Syndrome Virus (SFTS-NP) Antibody
<b>Catalog Number</b>	HM1156
<b>Description</b>	Ion Exchange Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Recombinant SFTS Virus nucleoprotein
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 90 % by SDS-PAGE
<b>Buffer</b>	0.15M Potassium Phosphate, 0.85% NaCl, pH 7.2
<b>Preservatives</b>	0.05% Sodium Azide