

Product Datasheet

Severe Fever with Thrombocytopenia Syndrome Virus (SFTS-NP) Antibody HM1156

Product Name	Severe Fever with Thrombocytopenia Syndrome Virus (SFTS-NP) Antibody
Catalog Number	HM1156
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant SFTS Virus nucleoprotein
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.15M Potassium Phosphate, 0.85% NaCl, pH 7.2
Preservatives	0.05% Sodium Azide