

# Product Datasheet

## Severe Fever with Thrombocytopenia Syndrome Virus (SFTS-NP) Antibody HM1157

Product Name	Severe Fever with Thrombocytopenia Syndrome Virus (SFTS-NP) Antibody
Catalog Number	HM1157
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant SFTS Virus nucleoprotein
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.15M Potassium Phosphate, 0.85% NaCl, pH 7.2
Preservatives	0.05% Sodium Azide